CONTENTS

English	2
Français	16
Deutsch	30
Español	44
Italiano	58

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➔ Visit our website : www.dentsplymallefer.com

Table of contents

Introduction.....	3
1. Indications for use	3
2. Contraindications	3
3. Warnings	3
4. Precautions	4
5. Adverse Reactions	4
6. Step by step instructions	5
A. Contents	5
B. Recharging the Battery	5
C. Replacement of the rechargeable battery	6
D. Getting Started	7
E. Search for the Apex	8
F. Sound Adjustment	10
G. Demo mode	10
H. Automatic Shutdown	11
I. Maintenance of your ProPex® II	11
J. Guarantee	11
K. Technical Specifications	12
L. Standard symbols	12
Disinfection, cleaning and sterilization	13

DIRECTIONS FOR USE

Introduction

Congratulations on the purchase of **ProPex® II**.

1. Indications for use

ProPex® II is an electronic device used for apex location determination and working length determination. This product must only be used in hospital environments, clinics or dental offices by qualified dental personnel.

2. Contraindications

The **ProPex® II** is not recommended for use in patients that have a pacemaker or other implanted electrical devices.

3. Warnings

See above.

- Only use the original charger

4. Precautions

The **ProPex® II** device must be used with the manufacturer's original accessories only.

In order to prevent infectious agent transmission it is highly recommended to use a rubber dam system during the endodontic procedure.

To ensure that short circuits do not impair the measurements, be particularly careful with patients fitted with metal crowns or bridges (avoid any metallic contact with the file or the lip clip)

- High concentrations of sodium hypochlorite may result in a lower accuracy. For working length determination, we recommend to use sodium hypochlorite solution at maximum 3% concentration.
- Make sure that the canal is wet enough to ensure reliability of the measurement
- Ensure that the file does not touch another instrument.
- Do not expose **ProPex® II** to any liquid.
- Avoid excessive liquids inside the tooth cavity to prevent overflow and incorrect measurements.
- **ProPex® II** must be stored in normal temperature (< 60 °C) and humidity conditions.
- Do not use in the presence of flammable products.
- Teeth with open apices will give imprecise results.
- If the bar graph suddenly makes a large movement in the upper part of the canal – continue slightly towards the apex so the signal returns normal.
- An apex locator may not be able to measure in all conditions. In any case, it is recommended to take an X-Ray prior to use of the unit and compare the information given by both means.
- For your own safety, please wear personal protective equipments (gloves, glasses, mask).

5. Adverse Reactions :

None known.

6. Step-by-Step Instructions

A. Content

Check the content of the equipment before use:

- one **ProPex[®] II** apex locator
- one charger
- one measurement cable
- two lip clips
- two connection hooks
- two connection forks
- one user manual

B. Recharging the Battery

The **ProPex[®] II** is delivered with a rechargeable battery.

On **ProPex[®] II** a screen indicator shows the level of the battery charge. When this is flashing, the battery requires recharging. However, it is still functional for several treatments before the battery shuts down.

Procedure for recharging the battery:

1. Disconnect the measurement cable
2. Connect the charger cable to the **ProPex[®] II**
3. Connect the charger to the mains

Charging screen will appear during battery charging. When the battery image on the screen stops blinking and turns green, the charging is complete.

Low

Mid

Full

Duration of charging: About 12 hours (24 hours after long periods without use).

Note: **ProPex[®] II** cannot be used while charging.

C. *Replacement of the rechargeable battery*

The battery compartment is located at the bottom of **ProPex® II** and its cover is secured by a screw.

1. Release the screw and remove the battery compartment cover
2. Remove the battery and disconnect the battery cable from the jack connector
3. Insert the new battery cable into the jack connector
4. Insert the battery into the battery compartment.
5. Close the battery compartment and secure it with the screw


Warning: Use only original battery pack from your supplier !

D. Getting Started

1. Disconnect the charger from the device if connected
2. Connect the measurement cable and switch the device On – after a short logo presentation, the main screen is displayed and the cursor on the tooth image will start blinking.
(While turning the device On without the measurement cable – the main screen will be displayed without cursor on the tooth image)
3. Attach the lip clip to the patient
4. Insert the file into the canal – (to ensure precise measurements the file size should be adjusted to the canal diameter)
5. Connect the hook to the file

The cursor on the tooth will stop blinking (accompanied by a double beep signal).


A blinking cursor and no sound signal indicates a faulty connection. Check your connections, clean your hook connection, moisten the canal, if necessary, and start again.

No other adjustments are necessary before starting measurement.

Note: Connection test feature is included in *ProPex® II* in order to check the cables:

- Connect metal part of the connection hook to the lip clip
- “Connection test” icon  should appear in the bottom left corner of the display, indicating proper connection
- If no icon appears, the connection hook or the measurement cable should be replaced


E. Searching for the Apex




1) Coronal and Medial Zone

Slowly introduce the file into the canal.

The cursor on the tooth icon indicates the progression of the file inside the canal.


Further progression of the file in the canal is shown both by the cursor on the tooth icon and by numerical value on the graphical scale (which does not represent length in millimeters). **ProPex® II** emits audible information of file progression via a series of progressive rate beeps.

2) Apical Zone

The apical zone is divided into 10 segments graduated from 0.9 to 0.0 (apex) as visual information of file progression.


When the apex is reached, the cursor shows “APEX” and a solid tone is emitted.

Once apex is reached, it is recommended to subtract 0.5 mm to determine the working length for shaping

3) Over-instrumentation

A red segment and a warning signal indicate that the file has passed the Apex. The cursor on the tooth icon shows “OVER”.

F. Sound Adjustment

ProPex® II is equipped with a sonic indicator which enables monitoring of the progression of the file within the canal.


This function, in conjunction with the display of progression, enables working "blind" while still monitoring the progression of the file.

The volume can be adjusted to one of four levels: mute, low, normal and high, by successive presses on the  volume key.

G. Demo mode

The built-in Demo mode is available to become acquainted with the device and to demonstrate it's operation.

1. Disconnect the measuring cable or the charger from the device if connected and turn Off the device
2. To start Demo Mode, press and hold the  (On/Off) button for about 2 sec. until the second beep sounds and "Demo" indication appears in the bottom left corner of the display

3. During Demo cycle device operating sequence is shown on the screen. Press  button to pause the simulation; press  button again to resume
4. When Demo cycle is completed, it is repeated automatically until interrupted by the operator
5. To exit Demo mode press  button and hold it for about 2 sec. until a beep sounds

Note: If measurement cable is connected to **ProPex® II** during Demo cycle, the device switches automatically to normal operation mode

H. Automatic Shutdown

ProPex[®] II automatically shuts down after 5 minutes without use. It is advisable, however, to manually switch off equipment after measurement by simply pressing the  key (On/Off).

I. Maintenance of your ProPex[®] II

The device does not contain user serviceable parts. The service and repair should be provided by factory trained service personnel only.

After each use, all the objects that were in contact with infectious agents should be cleaned using towels impregnated with a disinfecting and detergent solution (a bactericidal, fungicidal and aldehyde free solution). Use of chemical agents may cause damage to the equipment.

Reprocessing procedure :

- the lip clip, the hook and the fork must be sterilized between treatments by autoclaving at 134°C

See after for detailed procedure

The measuring cable can not be autoclaved.

J. Warranty

ProPex[®] II is warranted for 24 months from the date of purchase. The accessories (cables, battery etc.) are warranted for 6 months from the date of purchase.

The warranty is valid for normal usage conditions. Any modification or accidental damage will render the warranty void.

K. Technical Specifications

ProPex® II complies to IEC60601-1 safety standard and the requirement of CE Marking of Conformity.

ProPex® II electronic apex locator belongs to the following category of medical devices:

- Internally powered equipment (2.4V NiMH rechargeable battery)
- Type BF applied parts
- Not suitable for use in the presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide
- Continuous operation
- Ingress of liquids – not protected
- Environmental conditions during transportation: temperature: –20°C to +60°C (0 to 140°F); relative humidity: 10% to 90%, non-condensing

Technical specifications

Dimensions:	130 x 80 x 63 mm
Weight:	360 gr.
Type of screen:	Color Graphic TFT
Screen dimensions:	3.5 “
Supply:	2.4V NiMH rechargeable battery
External charger:	Input: 120 V / 50-60 Hz or 230V / 50-60 Hz Output: 6V DC

L. Standard symbols

On the device label appear standard symbols as follows:

Class II equipment

Type BF applied part

Attention, consult accompanying documents

Direct current

Recycling : PLEASE DO NOT THROW AWAY! This product and all its components must be recycled through your supplier.

Disinfection, cleaning and sterilization procedure **for the *ProPex*[®] II accessories**

Foreword

The lip clip, the hook and the fork must be cleaned, disinfected and sterilized before each use to prevent any contamination according to ISO 17664. This concerns the first use as well as the subsequent uses.

General recommendation

- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments.
- For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- Use only disinfecting solution which is approved for their efficacy (DGHM-listing, CE marking, FDA approval)

Procedure for lip clip, hook and fork

	Operation	Operating mode	Warning	
1.	Pre-Disinfection or Decontamination	- Soak immediately just after usage all instruments in a disinfectant solution combined with proteolytic enzyme if possible.	- Follow instructions and observe concentrations and immersion time given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments). - The disinfectant solution should be aldehyde free (to avoid blood impurities fixation). - Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (See general recommendation). - For visible impurities that are observed on instruments, a pre-cleaning is recommended by brushing them manually with soft material.	X
2.	Rinsing	- Rinse manually and abundantly the accessories with current water		X
3.	Manual Cleaning	- clean manually the accessories with an adequate brush, preliminary soaked in a clean pre-disinfectant solution	- the hook mechanism has to be activated during the cleaning process (press several times the push button) - no visible impurities should be observed on the accessories	X
4.	Rinsing	See point # 2		X
5.	Disinfection	- Immerse the accessories in a disinfectant solution (bactericidal, virucidal, fungicidal, tuberculocidal and aldehyde free) according to the manufacturer recommendations	- Follow instructions and observe concentrations and time given by the manufacturer	X
6.	Final rinsing	- See point # 2 - After rinsing, the accessories have to be dried.		X
7.	Inspection	- Inspect devices and sort out those with defects.	- Dirty instruments must be cleaned and disinfected again.	X
8.	Packaging	- Pack the devices in "Sterilization pouches".	- Check the validity period of the pouch given by the manufacturer to	X

			<p>determine the shelf life.</p> <ul style="list-style-type: none"> - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607. 	
9.	Sterilization	<ul style="list-style-type: none"> - Steam sterilization at: 134°C (237°F) during 18 min 	<ul style="list-style-type: none"> - The accessories (lip clip, hook and fork) must be sterilized according to the packaging labeling. - Use fractionated vacuum or gravity (less preferred) autoclaves (according to EN 13060, EN 285). - Use validated sterilization procedure according to ISO 17665-1 - Respect maintenance procedure of the autoclave device given by the manufacturer. - Use only the listed sterilization procedures. 	X
10.	Storage	<ul style="list-style-type: none"> - Keep devices in sterilization packaging in a dry and clean environment 	<ul style="list-style-type: none"> - Sterility cannot be guaranteed if packaging is open, damaged or wet (check the packaging before using the instruments). 	X

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