ORBI-Bond Universal Universal-Adhäsiv


EN **DE** **SV**

DA **NO** **FI**

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MP Medical Product GmbH
Graweg 18–22 | 27607 Geestland | Germany

REF 301967 ORBI-Bond Universal
Last revised: 2022-10
REF 301968 ORBI-Bond Universal Dual Cure Activator

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EN MD EU Medical device Instructions for use

Product description:

ORBI-Bond Universal is a dental adhesive for use with methacrylate-based restorative, core build-up or luting materials. It can be applied with phosphoric acid etching (selective enamel etching or total-etch) of the dental hard tissue as well as without (self-etch).

Indications:

- [A1] Adhesive for direct restorations of all classes of cavities using light-curing composite / compomer / Omrocer materials with a methacrylate base.
[A2] Adhesive for core build-ups made of light-curing composite materials.
[A3] Adhesive for indirect restorations with light-curing luting composites.
[B] Intraoral repair of composite, compomer or Omrocer restorations, ceramic veneers and all-ceramic restorations without an additional primer.

- [C] Treatment of hypersensitive tooth necks.
[D] Protective varnish for glass ionomer cement restorations.
[E1] Sealing of cavities prior to amalgam restorations.
[E2] Sealing of cavities and core preparations prior to the temporary luting of indirect restorations.

In combination with **ORBI-Bond Universal DC Activator (Dual Cure Activator)**:
[F1] Adhesive for direct self-curing or dual-curing composite restorations of all classes of cavities and for core build-ups.

- [F2] Adhesive for indirect restorations using dual-curing or self-curing luting composites.
[G] Luting of root posts with dual-curing or self-curing luting composites.

Contraindications:

ORBI-Bond Universal contains organic acids, (di)methacrylates and BHT. **DC Activator** contains amines. **ORBI-Bond Universal / DC Activator** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

ORBI-Bond Universal is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User: **ORBI-Bond Universal** should only be used by a professionally trained dental practitioner.

Use:

ORBI-Bond Universal not only ensures a secure bond to natural dental hard tissue, but also to composites, metals and ceramics.

It is necessary to isolate the operating site from moisture. Any contamination of the cavity or material surface with blood or saliva must be prevented. The use of a rubber dam is recommended.

The subsequent stages of application depend on the indication (see sections [A] - [G]). When using light-curing materials, **ORBI-Bond Universal** can be applied without the addition of **ORBI-Bond Universal Dual Cure Activator** (see sections [A] - [E]). The **ORBI-Bond Universal Dual Cure Activator** needs to be added only when working with self-curing or dual-curing materials (see sections [F] - [G]).

[A] Direct or indirect restorations with light-curing restorative, core build-up or luting composites

- Prepare the cavity and core according to the principles of the adhesive restoration technique, and remove any plaque or tartar in the case of noncarious class V cavities. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
- Etching of the dental hard tissue (optional)

In self-etch mode, **ORBI-Bond Universal** achieves high levels of adhesive strength on enamel and dentine. Prior etching of the dental hard tissue may optionally be carried out. Only enamel which is unprepared must always be etched at the beginning of the bonding process.

- Etch the surfaces to be etched using a dental acid-etch agent based on phosphoric acid (see applicable instructions for use).
- Aspirate the acid-etch agent, rinse with water for approx. 20 s. Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance.
- Place one drop of **ORBI-Bond Universal** on a mixing palette.
- Apply the adhesive evenly to the surfaces of the cavity and **rub it in** for 20 s with a disposable applicator.
- Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.
- Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).
- Apply the restorative, core build-up or luting material in accordance with the respective instructions for use.

[B] Intraoral repair of fillings or restorations (made from composite, glass ceramic, zirconium oxide ceramic, aluminium oxide ceramic, metal (incl. amalgam))

- Thoroughly roughen the surface of the restoration to which the adhesive is to be applied and sandblast it, preferably with aluminium oxide < 50 µm. In the case of ceramics, remove the varnish 1 mm beyond the margin.
- Clean thoroughly with water spray and dry with moisture- and oil-free air.
- Continue with the application of **ORBI-Bond Universal** immediately afterwards (see [A] 3. - 6.). The use of an additional primer is not necessary.
- Carry out the repair using a light-curing composite according to the manufacturer's instructions for use.

[C] Treatment of hypersensitive tooth necks

- Clean the tooth neck in the usual way (e.g. with polishing paste or a powder jet).
- Clean thoroughly with water spray and dry with moisture- and oil-free air. Do not overdry. The objective is a surface with a silky matt appearance.
- Afterwards, continue with the application of **ORBI-Bond Universal** (see [A] 3. - 6.).
- The oxygen inhibition layer should be removed using a cloth or a foam pellet.

[D] Protective varnish for glass ionomer cement restorations

- If a matrix was used to facilitate shaping, remove the matrix first.
- Place one drop of **ORBI-Bond Universal** on a mixing palette.
- Apply **ORBI-Bond Universal** to the modelled filling with a disposable applicator.
- Dry off **ORBI-Bond Universal** with dry, oil-free air for at least 5 s in order to remove any solvents.
- Cure **ORBI-Bond Universal** on every surface of the restoration using a commercially available polymerisation device for 10 s (LED or light with an output of > 500 mW/cm²).
- Once the restoration has completely cured, grind the occlusion and seal the ground areas once more using **ORBI-Bond Universal** in order o prevent the restoration from drying out.
- The oxygen inhibition layer should be removed using a cloth or a foam pellet.

[E] Sealing of cavities prior to amalgam restorations or prior to temporary luting
Note: **Do not** seal cavities with **ORBI-Bond Universal** if temporary luting using a methacrylate-based material is scheduled afterwards.

- Prepare the cavity according to the further treatment planned. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
- Place one drop of **ORBI-Bond Universal** on a mixing palette.
- Apply **ORBI-Bond Universal** to the cavity surface to be treated with a disposable applicator and rub it in for 20 s.
- Dry off **ORBI-Bond Universal** with dry, oil-free air for at least 5 s in order to remove any solvents.
- Cure **ORBI-Bond Universal** using a commercially available polymerisation device for 10 s (LED or halogen light with an output of > 500 mW/cm²).
- Depending on the indication, continue to work with the desired material afterwards in accordance with the corresponding instructions for use.

[F] Direct or indirect restorations with self- or dual-curing restorative, core build-up or luting composites

In combination with **ORBI-Bond Universal DC Activator, ORBI-Bond Universal** creates a strong bond to dual-curing and self-curing methacrylate-based composites and luting materials.

- Prepare the cavity or core according to the principles of the adhesive restoration technique. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
- Etching of the dental hard tissue (optional)

In the self-etch mode, **ORBI-Bond Universal** achieves high levels of adhesive strength on enamel and dentine. Prior etching of the dental hard tissue may optionally be carried out. Only enamel which is unprepared must always be etched at the beginning of the bonding process (e.g. in the luting of Maryland bridges).

- Etch the surfaces to be etched using a dental acid-etch agent based on phosphoric acid (see applicable instructions for use).
- Aspirate the acid-etch agent, rinse with water for approx. 20 s. Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance.

- Mix one drop of **ORBI-Bond Universal** and one drop of **ORBI-Bond Universal DC Activator** in a mixing palette with a disposable applicator for 3 s.
- Apply the mixed adhesive evenly on the surfaces of the cavity and rub it in for 20 s with a disposable applicator.
- Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.

- Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).
- Depending on the indication, continue to work with the desired material in accordance with the corresponding instructions for use.

[G] Luting of root posts

- Thoroughly mix one drop of **ORBI-Bond Universal** and one drop of **ORBI-Bond Universal DC Activator** in a mixing palette with a suitable disposable applicator (e.g. **ORBI-Brush**) for 3 s.
- Use the applicator to apply the mixed adhesive into the prepared root canal and to the occlusal surface of the core and **rub it in** for 20 s.
- Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents. Remove any excess material by dabbing with an absorbent paper strip. **Do not** light-cure the adhesive layer.
- Introduce the root post, which has been prepared in accordance with the manufacturer's instructions, into the root canal using a suitable self- or dualcuring composite for post luting (see applicable instructions for use) in such a way that a slight excess of material is achieved.
- To fix the root post and cure the occlusally-applied **ORBI-Bond Universal**, the

material should be light-cured for 40 s.

Warnings, precautionary measures:

- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of ORBI-Bond Universal. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with ORBI-Bond Universal.
- Avoid contact with the oral mucosa. Contact with **ORBI-Bond Universal** may cause whitening of the oral mucosa as a result of the coagulation of proteins. This is a reversible effect which subsides within a few days.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- ORBI-Bond Universal** is light-curing. Hence, please avoid exposure to strong lighting during application.
- The mixed adhesive is self-curing and light-curing, further processing must therefore take place immediately after mixing.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
ORBI-Bond Universal: BisGMA, HEMA, ethanol, water, HEDMA, methacrylate phophoric acid ester, methacrylate functionalised polyacid, UDMA, initiators, stabilisers
ORBI-Bond Universal and DC Activator: Ethanol, BisGMA, HEMA, water, HEDMA, methacrylate phophoric acid ester, methacrylate functionalised polyacid, UDMA, initiators, stabilisers

Storage:
Tightly close **ORBI-Bond Universal** and **ORBI-Bond Universal DC Activator** immediately after use.

Store at temperatures of 2 °C to 8 °C (**ORBI-Bond Universal**) or direct 4 °C to 23 °C (**ORBI-Bond Universal Dual Cure Activator**). Do not expose to higher sunlight or the operating light. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **ORBI-Bond Universal / ORBI-Bond Universal DC Activator** must be reported to the manufacturer and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of **ORBI-Bond Universal / ORBI-Bond Universal DC Activator** can be found in the European database on medical devices (EUDAMED – https://ec.europa.eu/tools/eudamed).

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