

SAFETY DATA SHEET

(Revision:7/5/2022; Supersedes: 6/3/2022)

gulation (EC) No 1907/2006 (REACH)				
Section 1 Identification of th	e Substance/Prepara	tion and of the Company/U	ndertaking.	
1.1 Product Identifier Product Type: Model \$	Stones. Plasters and D	Die Materials		
Trade Names:				
Bitestone Jade Stone FlowStone, Black Prima-Rock Mounting Stone Ulti Rock PL Lab Plaster Fast Set	Buffstone Hard Rock Lean Rock Ivory Quickstone Silky-Rock EuroStone PL Lab Stone Blue	Die Stone, Ivory Die Stone HS Microstone RapidFlask Snap Stone Economy Stone PL DPG Labstone Buff	FlowStone Laboratory Plaster Mounting Plaster ResinRock Orthodontic Plaster* Orthodontic Stone* Super Die	
1.2 Relevant Identified Uses Product Use: Modeling st Uses Advised Against: F	ones and plasters		Against	
1.3 Details of the Supplier of	the Substance or Mi	xture		
Manufacturer: Whip Mix Corporation 361 Farmington Avenue Louisville, Kentucky, USA Emergency Telephone Num Fax Number: (502) 634-451	nber: (502) 634-1451	EU Importer Whip Mix Europe GmbH Wißstrasse 26 – 28 D – 44137 Dortmund, G +49 (0) 231 / 567 70 8-0 +49 (0) 231 / 567 70 850	ermany	
1.4 Emergency Telephone N Transportation Emergence	es: CHEMTREC 1(d	800) 424-9300 (U.S. and Canad Ils: 1- 703-527-3887 (Collect ca		
Other Product Information: www.whipmix.com Email: Info@whipmix.com				
Section 2 Hazard Identification	on.			
2.1 Classification of the Substa	ance or Mixture:			
OSHA/WHMIS/GHS/CLP (Classification (1272/2			
Health Hazards		Physical Hazards	Environmental Hazards	
Not Hazardous		Not Hazardous	Not Hazardous	
2.2 Label Elements None required				
2.3 Other Hazards: None				

Section 3 Composition/Information on Ingredients.

<u>Substance</u>	<u>CAS No. /</u> <u>EC Number</u>	<u>%</u>	CLP/GHS Classification (1272/2008)
Plaster of Paris (Calcium Sulfate Hemihydrate)	26499-65-0 / 607-950-0	90 – 100	Not hazardous

See Section 16 for full text of GHS Classifications.

Section 4 First-Aid Measures.

4.1 Description of First Aid Measures

Inhalation: Remove exposed person to fresh air. If irritation or other symptoms persist, get medical attention. **Eyes:** Flush with large quantities of water, holding the eyelids apart. If irritation persists consult a physician.

Skin: No first aid is generally required. Wash skin with soap and water.

Ingestion: May cause gastrointestinal discomfort and intestinal blockage. If swallowed, drink 1 or 2 glasses of water to dilute. Never give anything by mouth to an unconscious or convulsing person. Get immediate medical attention.

4.2 Most Important symptoms and effects, both acute and delayed: May cause eye irritation. Inhalation of dust may cause mucous membrane and respiratory irritation. When mixed with water, this material hardens and becomes very hot – may cause burns.

4.3 Indication of any immediate medical attention and special treatment needed: Immediate medical attention is required for ingestions.

Section 5 Fire-Fighting Measures.

5.1 Extinguishing Media: Use media appropriate for surrounding fire. Water may cause product to solidify.
5.2 Special Hazards Arising from the Substance or Mixture: The product does not burn but will decompose producing calcium oxide and sulfur oxides.

5.3 Advice for Fire-Fighters: Firefighters should wear full emergency equipment and approved positive pressure self-contained breathing apparatus. Cool fire exposed containers with water.

Section 6 Accidental Release Measures.

6.1 Personal Precautions, Protective Equipment and Emergency Procedures: Wear appropriate protective clothing as described in Section 8.

6.2 Environmental Precautions: Report releases as required by local and national authorities.

6.3 Methods and Material for Containment and Cleaning Up: Collect using dustless method (HEPA vacuum or wet method) and place in appropriate container for use. Do not use compressed air.

6.4 Reference to Other Sections: Refer to Section 8 for personal protective equipment and Section 13 for disposal information.

Section 7 Handling and Storage.

7.1 Precautions for Safe Handling: Avoid contact with eyes. Do not breathe dust. Wear protective clothing and equipment as described in Section 8. Use with adequate ventilation and proper dust collection methods to keep exposure level below occupational exposure limits. Wash thoroughly with soap and water after handling. Keep containers closed when not in use.

7.2 Conditions for Safe Storage, Including any Incompatibilities: Store in a cool, dry, well-ventilated area away from incompatible materials. Protect from physical damage.

7.3 Specific end use(s):

Industrial uses: None identified

Professional uses: Model stones, plaster and die materials for dental technicians.

Section 8 Exposure Controls/Personal Protection.

8.1 Control Parameters:

Plaster of Paris (Calcium Sulfate	15 mg/m ³ (total dust); 5 mg/m ³ (respirable fraction)
Hemihydrate)	OSHA PEL

Refer to local regulations for exposure limits not listed above

8.2 Exposure Controls: Recommended Monitoring Procedures: None.

Appropriate engineering controls: Use with adequate local exhaust ventilation to maintain exposures below the occupational exposure limits.

Personal Protective Measurers

Respiratory protection: If the exposure levels are exceeded and irritation is experienced an approved dust/mist respirator appropriate for the form and concentration of the contaminants should be used. In the USA refer to OSHA regulations, in the EU refer to EN Standards (EN 149 or 405). Selection and use of respiratory equipment must be in accordance with applicable regulations and good industrial hygiene practice.

Skin protection: For prolonged use or in dusty conditions, wear rubber gloves. In the EU refer to EN 374. **Eye protection:** Chemical safety goggles if needed to avoid eye contact. In the EU refer to EN 166. **Other:** Impervious clothing as needed to avoid contamination of personal clothing.

Section 9 Physical and Chemical Properties.

9.1 Information on basic Physical and Chemical Properties

Appearance: Powder, with variety of colors Physical State: Solid Odor: Odorless.

Odor threshold: Not applicable Melting point/freezing point: 293°F / 145°C Flash point: Not applicable Flammability (solid, liquid, gas): Not applicable Flammable limits: LEL: Not applicable Vapor pressure: Not applicable Relative density: 2.5 – 3.5 Partition coefficient: n-octanol/water: Not available Decomposition temperature: 2642°F / 1450°C Explosive Properties: Not applicable Particle Characteristics: No data available pH: Not availableBoiling point: Not applicableEvaporation rate: Not applicable

UEL: Not applicable Relative vapor density (air = 1): Not applicable Solubility In Water: 0.2% Auto-ignition temperature: Not applicable

Kinematic Viscosity: Not applicable Oxidizing Properties: Not applicable

9.2 Other Information: None available

Section 10 Stability and Reactivity.

10.1 Reactivity: None known.

10.2 Chemical stability: Stable

10.3 Possibility of hazardous reactions: None known.

10.4 Conditions to avoid: Avoid unintentional contact with water. Product will harden and produce heat.

10.5 Incompatible materials: Avoid acids and oxidizing agents.

10.6 Hazardous decomposition products: Thermal decomposition (above 2642°F/1450°C) may generate calcium oxide and sulfur dioxide.

Section 11 Toxicological Information.

11.1 Information on Toxicological Effects:

Potential Health Effects:

Eyes: Dust may cause mechanical irritation and possible injury. **Skin:** Dust may cause irritation. When mixed with water, the plaster of paris hardens and becomes hot – may cause skin burns. **Ingestion:** No adverse effects expected for normal, incidental ingestion. Large amounts may cause gastrointestinal blockage and discomfort.

Inhalation: Inhalation of dust may cause irritation to the nose, throat and upper respiratory tract with coughing and shortness of breath.

Chronic Health Effects: None known.

Acute Toxicity: Based on available data, the classification criteria are not met. Plaster of Paris: Oral rat LD50 > 2000 mg/kg; Inhalation rat LC50 > 3.26 mg/L/4 hr (structurally similar chemical)

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Eye damage/ irritation: Based on available data, the classification criteria are not met.

Respiratory Irritation: Based on available data, the classification criteria are not met.

Respiratory Sensitization: Based on available data, the classification criteria are not met.

Skin Sensitization: Based on available data, the classification criteria are not met.

Germ Cell Mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met. None of the components of this product are listed as carcinogens by OSHA, IARC or NTP or the EU CLP.

Reproductive Toxicity: Based on available data, the classification criteria are not met. **Specific Target Organ Toxicity**:

Single Exposure: Based on available data, the classification criteria are not met.. Repeat Exposure: Based on available data, the classification criteria are not met.

Aspiration Toxicity: Based on available data, the classification criteria are not met.

11.2 Information on Other Hazards

11.2.1 Endocrine Disrupting Properties: None known

11.2.2 Other Information: None

Section 12 Ecological Data.

12.1 Ecotoxicity:

Plaster of Paris: 96 hr LC50 Pimephales promelas>1970 mg/L, 48 hr LD50 daphnia magna >79 mg/L, 72 hr EC50 Pseudokirchnerella subcapitata >79 mg/L (structurally similar chemical)

12.2 Persistence and degradability: Not applicable

12.3 Bioaccumulative potential: No data available

12.4 Mobility in soil: No data available

12.5 Results of PBT and vPvB assessment: Components do not meet the criteria of PBT or vPvB.

12.6 Endocrine Disrupting Properties: None known

12.7 Other adverse effects: Not required.

Section 13 Disposal Considerations.

13.1 Waste Treatment Methods: Dispose in accordance with all national and local regulations.

Section 14. Transport Information.

	14.1 UN Number	14.2 UN Proper Shipping Name	14.3 Hazard Class(s)	14.4 Packing Group	14.5 Environmental Hazards
US DOT		Not Regulated			
Canadian TDG		Not Regulated			
EU ADR/RID		Not Regulated			
IMDG		Not Regulated			
IATA/ICAO		Not Regulated			

14.6 Special precautions for User: Not applicable

14.7 Transport in Bulk According to IMO Instruments: Not applicable – product is transported only in packaged form.

Section 15 Regulatory Information.

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture

US Regulations

SARA Section 313 (40 CFR 372): This product contains the following toxic chemical(s) subject to reporting requirements of SARA 313: None

SARA Section 311/312 (40 CFR 370) Hazard Categories: Refer to Section 2 for the OSHA Hazard Classification.

Comprehensive Environmental Response and Liability Act of 1980 (CERCLA): This product is not subject to CERCLA reporting requirements. Many states have more stringent release reporting requirements. Report spills required under federal, state and local regulations.

Toxic Substances Control Act (TSCA): The components of this product are not listed on the TSCA inventory

California: This product contains the following substances known to the state of California to cause cancer and/or reproductive toxicity: None known

Canadian Regulations

Canadian Environmental Protection Act: All of the components of this product are listed on the Canadian Domestic Substances List (DSL) or has a quantity limitation.

International Chemical Inventories

Australia: Information not available for this product.

China: All of the components in this product are listed on the Inventory of Existing Chemical Substances in China (IECSC) or exempt.

European Union: All the components in this product are listed on the EINECS inventory or exempt.

Korea: Information not available for this product.

New Zealand: One or more components in this product are listed on the New Zealand Inventory of Chemicals (NZIoC) or exempt.

Philippines: Information not available for this product.

15.2 Chemical safety assessment: None required.

16. Other Information.

HMIS Rating: Health 1 Flammability 0 Physical Hazard 0 Hazard: 4-Severe; 3-Serious; 2-Moderate; 1-Slight; 0-Minimum

<u>CLP/GHS Classification and H Phrases for Reference (See Section 3)</u> None.

Prepared By: Denuse A. Deids	Translated By:
Date: July 5, 2022	Date: