




# SpeedCEM® Plus

Placement of crowns and bridges (oxide ceramic, metal-ceramic, lithium-disilicate)  
Befestigung von Kronen und Brücken (Oxidkeramik, Metallkeramik, Lithium-Disilikat)


Zirconium oxide, aluminium oxide, metal-ceramic, metal.  
Zirkoniumoxid, Aluminiumoxid, Metallkeramik, Metall.




Sandblast  
Sandstrahlen




Ivoclean®







20s







Lithium-disilicate.  
Lithium-Disilikat.



Monobond  
Etch & Prime®



10s



CE 0123

Rx only



Manufacturer:  
Ivoclar Vivadent AG, 9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)



see instructions

ivoclar  
vivadent  
clinical